

SLMAP antibody

Product Information

Catalog No.:	FNab07986
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role during myoblast fusion.

Immunogen information

Immunogen:	sarcolemma associated protein
Synonyms:	Sarcolemmal membrane-associated protein (Sarcolemmal-associated protein) SLMAP KIAA1601 SLAP
Observed MW:	35, 45, 63-66, 83-91 kDa
Uniprot ID :	Q14BN4

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab07986(SLMAP Antibody) at dilution of 1:600