

SMAP1 antibody

Product Information

Catalog No.:	FNab08002
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


GTPase activating protein that acts on ARF6. Plays a role in clathrin-dependent endocytosis. May play a role in erythropoiesis(By similarity).

Immunogen information


Immunogen:	small ArfGAP 1
Synonyms:	Stromal membrane-associated protein 1 SMAP1
Observed MW:	50-55 kDa
Uniprot ID :	Q8IYB5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach tissue slide using FNab08002(SMAP1 Antibody) at dilution of 1:50



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab08002(SMAP1 antibody) at dilution of 1:1000