

SMARCA4 antibody

Product Information


Catalog No.:	FNab08005
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 4
Synonyms:	Transcription activator BRG1 ATP-dependent helicase SMARCA4 BRG1-associated factor 190A (BAF190A) Mitotic growth and transcription activator Protein BRG-1 Protein brahma homolog 1 SNF2-beta SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily A member 4 SMARCA4 BAF190A BRG1 SNF2B SNF2L4
Observed MW:	185 kDa
Uniprot ID :	P51532

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP:1:500-1:2000; IHC: 1:50-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab08005(SMARCA4 antibody) at dilution of 1:100



IP Result of anti-SMARCA4 (IP: FNab08005, 5ug; Detection: FNab08005 1:1000) with HeLa cells lysate 2500ug



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08005(SMARCA4 antibody) at dilution of 1:800