

SMARCD2 antibody

Product Information

Catalog No.:	FNab08012
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in transcriptional activation and repression of select genes by chromatin remodeling(alteration of DNA-nucleosome topology).

Immunogen information


Immunogen:	SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily d, member 2
Synonyms:	SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily D member 2 60 kDa BRG-1/Brm-associated factor subunit B BRG1-associated factor 60B (BAF60B) SMARCD2 BAF60B
Observed MW:	52 kDa
Uniprot ID :	Q92925

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab08012(SMARCD2 antibody) at dilution of 1:50



RAW264.7 cells were subjected to SDS PAGE followed by western blot with FNab08012(SMARCD2 antibody) at dilution of 1:500