

## PPP4R3B antibody

### Product Information

Catalog No.:	FNab08024
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	SMEK homolog 2, suppressor of mek1(Dictyostelium)
Synonyms:	Serine/threonine-protein phosphatase 4 regulatory subunit 3B SMEK homolog 2 PPP4R3B KIAA1387 PP4R3B SMEK2
Observed MW:	51 kDa
Uniprot ID :	Q5MIZ7

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab08024(SMEK2 antibody) at dilution of 1:100



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab08024(SMEK2 antibody) at dilution of 1:500