

## MYH11 antibody

### Product Information

Catalog No.:	FNab08026
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	myosin, heavy chain 11, smooth muscle
Synonyms:	Myosin-11 Myosin heavy chain 11 Myosin heavy chain, smooth muscle isoform SMMHC MYH11 KIAA0866
Observed MW:	230 kDa
Uniprot ID :	P35749

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse uterus tissue were subjected to SDS PAGE followed by western blot with FNab08026(SMMHC antibody) at dilution of 1:3000