

MYH11 antibody

Product Information

Catalog No.:	FNab08027
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	myosin, heavy chain 11, smooth muscle
Synonyms:	Myosin-11 Myosin heavy chain 11 Myosin heavy chain, smooth muscle isoform SMMHC MYH11 KIAA0866
Observed MW:	230 kDa
Uniprot ID :	P35749

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle using FNab08027(SMMHC antibody) at dilution of 1:50



A375 cells were subjected to SDS PAGE followed by western blot with FNab08027(SMMHC antibody) at dilution of 1:1000