

SMPX antibody

Product Information

Catalog No.:	FNab08042
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in the regulatory network through which muscle cells coordinate their structural and functional states during growth, adaptation, and repair.

Immunogen information


Immunogen:	small muscle protein, X-linked
Synonyms:	Small muscular protein Stretch-responsive skeletal muscle protein SMPX SRMX
Observed MW:	9 kDa
Uniprot ID :	Q9UHP9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab08042(SMPX antibody) at dilution of 1:50



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab08042(SMPX antibody) at dilution of 1:400