

SMU1 antibody

Product Information

Catalog No.:	FNab08045
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	smu-1 suppressor of mec-8 and unc-52 homolog(C. elegans)
Synonyms:	WD40 repeat-containing protein SMU1 Smu-1 suppressor of mec-8 and unc-52 protein homolog WD40 repeat-containing protein SMU1, N-terminally processed SMU1
Observed MW:	45-50 kDa
Uniprot ID :	Q2TAY7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab08045(SMU1 antibody) at dilution of 1:100



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab08045(SMU1 antibody) at dilution of 1:1000