

SNAI3 antibody

Product Information

Catalog No.:	FNab08053
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Seems to inhibit myoblast differentiation. Transcriptional repressor of E-box-dependent transactivation of downstream myogenic bHLHs genes. Binds preferentially to the canonical E-box sequences 5'-CAGGTG-3' and 5'-CACCTG-3'(By similarity).

Immunogen information

Immunogen:	snail homolog 3(Drosophila)
Synonyms:	Zinc finger protein SNAI3 Protein snail homolog 3 Zinc finger protein 293 SNAI3 ZNF293
Observed MW:	32-38 kDa
Uniprot ID :	Q3KNW1

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab08053(SNAI3 antibody) at dilution of 1:1000