

SNRPA1 antibody

Product Information

Catalog No.:	FNab08069
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


This protein is associated with sn-RNP U2. It helps the A' protein to bind stem loop IV of U2 snRNA.

Immunogen information


Immunogen:	small nuclear ribonucleoprotein polypeptide A'
Synonyms:	U2 small nuclear ribonucleoprotein A' (U2 snRNP A') SNRPA1
Observed MW:	28 kDa
Uniprot ID :	P09661

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung using FNab08069(SNRPA1 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08069(SNRPA1 antibody) at dilution of 1:500