

SNRPD2 antibody

Product Information

Catalog No.:	FNab08073
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


The protein encoded by this gene belongs to the small nuclear ribonucleoprotein core protein family. It is required for pre-mRNA splicing and small nuclear ribonucleoprotein biogenesis. Multiple transcript variants encoding different isoforms have been found for this gene.

Immunogen information


Immunogen:	small nuclear ribonucleoprotein D2 polypeptide 16.5kDa
Synonyms:	Small nuclear ribonucleoprotein Sm D2 (Sm-D2) snRNP core protein D2 SNRPD2 SNRPD1
Observed MW:	13 kDa
Uniprot ID :	P62316

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNa08073(SNRNP2 Antibody) at dilution of 1:50



A2780 cells were subjected to SDS PAGE followed by western blot with FNa08073(SNRNP2 antibody) at dilution of 1:1000