

SNRPN antibody

Product Information

Catalog No.:	FNab08079
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in tissue-specific alternative RNA processing events.

Immunogen information


Immunogen:	small nuclear ribonucleoprotein polypeptide N
Synonyms:	Small nuclear ribonucleoprotein-associated protein N (snRNP-N) Sm protein D (Sm-D) Sm protein N (Sm-N, SmN) Tissue-specific-splicing protein SNRPN HCERN3 SMN
Observed MW:	25-30 kDa
Uniprot ID :	P63162

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab08079(SNRPN antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab08079(SNRPN antibody) at dilution of 1:800