

SNX1 antibody

Product Information

Catalog No.:	FNab08082
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This endosomal protein regulates the cell-surface expression of epidermal growth factor receptor. This protein also has a role in sorting protease-activated receptor-1 from early endosomes to lysosomes. This protein may form oligomeric complexes with family members. This gene results in three transcript variants encoding distinct isoforms.

Immunogen information


Immunogen:	sorting nexin 1
Synonyms:	Sorting nexin-1 SNX1
Observed MW:	80 kDa
Uniprot ID :	Q13596

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000 - 1:2000; IHC: 1:100 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab08082(SNX1 Antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08082(SNX1 antibody) at dilution of 1:1000