

SNX4 antibody

Product Information

Catalog No.:	FNab08093
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in several stages of intracellular trafficking. Plays a role in recycling endocytosed transferrin receptor and prevent its degradation.

Immunogen information

Immunogen:	sorting nexin 4
Synonyms:	Sorting nexin-4 SNX4
Observed MW:	50 kDa
Uniprot ID :	O95219

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab08093(SNX4 antibody) at dilution of 1:300