

## SST antibody

### Product Information

Catalog No.:	FNab08109
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	somatostatin
Synonyms:	Somatostatin Growth hormone release-inhibiting factor Somatostatin-28 Somatostatin-14 (SST-14) Neuronostatin (NST) SST
Observed MW:	17 kDa
Uniprot ID :	P61278

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:500
Image:	



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab08109(SST Antibody) at dilution of 1:200



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab08109(SST antibody) at dilution of 1:1000