

## SORT1 antibody

### Product Information

Catalog No.:	FNab08117
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	sortilin 1
Synonyms:	Sortilin 100 kDa NT receptor Glycoprotein 95 (Gp95) Neurotensin receptor 3 (NT3, NTR3) SORT1
Observed MW:	100-110 kDa
Uniprot ID :	Q99523

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab08117(SORT1 Antibody) at dilution of 1:100



IP Result of anti-Sortilin (IP:FNab08117, 4ug; Detection:FNab08117 1:500) with rat brain tissue lysate 4000ug.



human brain tissue were subjected to SDS PAGE followed by western blot with FNab08117(Sortilin antibody) at dilution of 1:1200