

SOX15 antibody

Product Information

Catalog No.:	FNab08123
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds to the 5'-AACAAT-3' sequence.

Immunogen information

Immunogen:	SRY(sex determining region Y)-box 15
Synonyms:	Protein SOX-15 Protein SOX-12 Protein SOX-20 SOX15 SOX12 SOX20 SOX26 SOX27
Observed MW:	30 kDa
Uniprot ID :	O60248

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab08123(SOX15 Antibody) at dilution of 1:300