

## SOX6 antibody

### Product Information

Catalog No.:	FNab08128
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Transcriptional activator. Binds specifically to the DNA sequence 5'-AACCAAT-3'. Plays a key role in several developmental processes, including neurogenesis and skeleton formation.

### Immunogen information

Immunogen:	SRY(sex determining region Y)-box 6
Synonyms:	Transcription factor SOX-6 SOX6
Observed MW:	90 kDa
Uniprot ID :	P35712

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab08128(SOX6 antibody) at dilution of 1:500