

## SOX8 antibody

### Product Information

Catalog No.:	FNab08130
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in central nervous system, limb and facial development. May be involved in male sex determination. Binds the consensus motif 5'-[AT][AT]CAA[AT]G-3'(By similarity).

### Immunogen information

Immunogen:	SRY(sex determining region Y)-box 8
Synonyms:	Transcription factor SOX-8 SOX8
Observed MW:	48 kDa
Uniprot ID :	P57073

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab08130(SOX8 antibody) at dilution of 1:500