

SPAG11A antibody

Product Information

Catalog No.:	FNab08139
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	sperm associated antigen 11A
Synonyms:	Sperm-associated antigen 11A Antimicrobial-like protein Bin-1b Human epididymis-specific protein 2 (He2) Protein EP2 Sperm antigen HE2 SPAG11A BIN1B EP2 HE2
Observed MW:	20 kDa
Uniprot ID :	Q6PDA7

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab08139(SPAG11A antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab08139(SPAG11A antibody) at dilution of 1:300