

SUN5 antibody

Product Information

Catalog No.:	FNab08141
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	sperm associated antigen 4-like
Synonyms:	SUN domain-containing protein 5 Sad1 and UNC84 domain-containing protein 5 Sperm-associated antigen 4-like protein Testis and spermatogenesis-related gene 4 protein SUN5 SPAG4L TSARG4
Observed MW:	43 kDa
Uniprot ID :	Q8TC36

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



human testis tissue were subjected to SDS PAGE followed by western blot with FNab08141(SPAG4L antibody) at dilution of 1:400