

## SPAG4-Specific antibody

### Product Information

Catalog No.:	FNab08142
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	sperm associated antigen 4
Synonyms:	Sperm-associated antigen 4 protein Outer dense fiber-associated protein SPAG4 SUN domain-containing protein 4 SPAG4 SUN4
Observed MW:	55 kDa
Uniprot ID :	Q9NPE6

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse testis tissue slide using FNab08142(SPAG4-Specific Antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab08142(SPAG4-Specific Antibody) at dilution of 1:600