

SPANXA1 antibody

Product Information

Catalog No.:	FNab08149
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	SPANX family, member A2
Synonyms:	Sperm protein associated with the nucleus on the X chromosome A Cancer/testis antigen 11.1 (CT11.1) Nuclear-associated protein SPAN-Xa (SPAN-X, SPANX-A) SPANX family member A SPANXA1 SPANXA
Observed MW:	27 kDa
Uniprot ID :	Q9NS26

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab08149 (SPANXA2 Antibody) at dilution of 1:50



DU 145 cells were subjected to SDS PAGE followed by western blot with FNab08149 (SPANXA2 antibody) at dilution of 1:500