

SPATA3 antibody

Product Information

Catalog No.:	FNab08163
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	spermatogenesis associated 3
Synonyms:	Spermatogenesis-associated protein 3 Testis and spermatogenesis cell-related protein 1 Testis spermatocyte apoptosis-related protein 1 SPATA3 TSARG1
Observed MW:	None
Uniprot ID :	Q8NHX4

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab08163 (SPATA3 Antibody) at dilution of 1:50