

SPATA7 antibody

Product Information

Catalog No.:	FNab08165
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in photoreceptor cells maintenance(By similarity). It is required for recruitment and proper localization of RPGRIP1 to the photoreceptor connecting cilium(CC), as well as protein trafficking across the CC to the outer segments(By similarity).

Immunogen information


Immunogen:	spermatogenesis associated 7
Synonyms:	Spermatogenesis-associated protein 7 HSD-3.1 Spermatogenesis-associated protein HSD3 SPATA7 HSD3
Observed MW:	68-80 kDa
Uniprot ID :	Q9P0W8

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IP, IHC
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab08165(SPATA7 Antibody) at dilution of 1:50



IP Result of anti-SPATA7 (IP:FNab08165, 3ug; Detection:FNab08165 1:300) with mouse testis tissue lysate 8000ug.



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08165(SPATA7 antibody) at dilution of 1:300