

SPATS2L antibody

Product Information

Catalog No.:	FNab08166
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	viral DNA polymerase-transactivated protein 6
Synonyms:	SPATS2-like protein DNA polymerase-transactivated protein 6 Stress granule and nucleolar protein (SGNP) SPATS2L DNAPTP6
Observed MW:	54 kDa
Uniprot ID :	Q9NUQ6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab08166(SPATS2L antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab08166(SPATS2L antibody) at dilution of 1:600