

## SPIN2B antibody

### Product Information

Catalog No.:	FNab08179
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


### Background

#### Immunogen information


Immunogen:	spindlin family, member 2B
Synonyms:	Spindlin-2B Spindlin-like protein 2B (SPIN-2, SPIN-2B) SPIN2B SPIN2
Observed MW:	29 kDa
Uniprot ID :	Q9BPZ2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab08179(SPIN2B antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab08179(SPIN2B antibody) at dilution of 1:1000