

SPINK1 antibody

Product Information

Catalog No.:	FNab08180
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Serine protease inhibitor which exhibits anti-trypsin activity(PubMed:7142173). In the pancreas, protects against trypsin-catalyzed premature activation of zymogens(By similarity). In the male reproductive tract, binds to sperm heads where it modulates sperm capacitance by inhibiting calcium uptake and nitrogen oxide(NO) production.

Immunogen information

Immunogen:	serine peptidase inhibitor, Kazal type 1
Synonyms:	Serine protease inhibitor Kazal-type 1 Pancreatic secretory trypsin inhibitor Tumor-associated trypsin inhibitor (TATI) SPINK1 PSTI
Observed MW:	None
Uniprot ID :	P00995

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab08180(SPINK1 antibody) at dilution of 1:200