

SPOCK2 antibody

Product Information

Catalog No.:	FNab08187
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May participate in diverse steps of neurogenesis. Binds calcium.

Immunogen information

Immunogen:	sparc/osteonectin, cwcv and kazal-like domains proteoglycan(testican) 2
Synonyms:	Testican-2 SPARC/osteonectin, CWCV, and Kazal-like domains proteoglycan 2 SPOCK2 KIAA0275 TICN2
Observed MW:	60 kDa
Uniprot ID :	Q92563

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08187 (SPOCK2 antibody) at dilution of 1:1000