

## SPON2 antibody

### Product Information

Catalog No.:	FNab08188
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Cell adhesion protein that promotes adhesion and outgrowth of hippocampal embryonic neurons. Binds directly to bacteria and their components and functions as an opsonin for macrophage phagocytosis of bacteria. Essential in the initiation of the innate immune response and represents a unique pattern-recognition molecule in the ECM for microbial pathogens(By similarity). Binds bacterial lipopolysaccharide(LPS).

### Immunogen information

Immunogen:	spondin 2, extracellular matrix protein
Synonyms:	Spondin-2 Differentially expressed in cancerous and non-cancerous lung cells 1 (DIL-1) Mindin SPON2 DIL1
Observed MW:	None
Uniprot ID :	Q9BUD6

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using FNab08188(SPON2 antibody) at dilution of 1:50