

HM13 antibody

Product Information

Catalog No.:	FNab08190
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	histocompatibility(minor) 13
Synonyms:	Minor histocompatibility antigen H13 Intramembrane protease 1 (IMP-1, IMPAS-1, hIMP1) Presenilin-like protein 3 Signal peptide peptidase HM13 H13 IMP1 PSL3 SPP
Observed MW:	45 kDa, 100 kDa
Uniprot ID :	Q8TCT9

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab08190(SPP antibody) at dilution of 1:100



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab08190(SPP antibody) at dilution of 1:300