

SPTBN1 antibody

Product Information

Catalog No.:	FNab08203
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Fodrin, which seems to be involved in secretion, interacts with calmodulin in a calcium-dependent manner and is thus candidate for the calcium-dependent movement of the cytoskeleton at the membrane.

Immunogen information

Immunogen:	spectrin, beta, non-erythrocytic 1
Synonyms:	Spectrin beta chain, non-erythrocytic 1 Beta-II spectrin Fodrin beta chain Spectrin, non-erythroid beta chain 1 SPTBN1 SPTB2
Observed MW:	275 kDa
Uniprot ID :	Q01082

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab08203(SPTBN1 Antibody) at dilution of 1:600