

## SPTBN2 antibody

### Product Information

Catalog No.:	FNab08205
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Probably plays an important role in neuronal membrane skeleton.

### Immunogen information

Immunogen:	spectrin, beta, non-erythrocytic 2
Synonyms:	Spectrin beta chain, non-erythrocytic 2 Beta-III spectrin Spinocerebellar ataxia 5 protein SPTBN2 KIAA0302 SCA5
Observed MW:	275 kDa
Uniprot ID :	O15020

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab08205(SPTBN2 antibody) at dilution of 1:500