

SCARB1 antibody

Product Information

Catalog No.:	FNab08215
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	scavenger receptor class B, member 1
Synonyms:	Scavenger receptor class B member 1 (SRB1) CD36 and LIMPII analogous 1 (CLA-1) CD36 antigen-like 1 Collagen type I receptor, thrombospondin receptor-like 1 SR-BI SCARB1 CD36L1 CLA1
Observed MW:	82 kDa
Uniprot ID :	Q8WTV0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:1000; IF: 1:10-1:100
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab08215(SR-BI antibody) at dilution of 1:500