

NCOA1 antibody

Product Information

Catalog No.:	FNab08217
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	nuclear receptor coactivator 1
Synonyms:	Nuclear receptor coactivator 1 (NCoA-1) Class E basic helix-loop-helix protein 74 (bHLHe74) Protein Hin-2 RIP160 Renal carcinoma antigen NY-REN-52 Steroid receptor coactivator 1 (SRC-1) NCOA1 BHLHE74 SRC1
Observed MW:	160 kDa
Uniprot ID :	Q15788

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab08217(NCOA1 antibody) at dilution of 1:100



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08217(NCOA1 antibody) at dilution of 1:400