

NCOA3 antibody

Product Information

Catalog No.:	FNab08218
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	nuclear receptor coactivator 3
Synonyms:	Nuclear receptor coactivator 3 (NCoA-3) ACTR Amplified in breast cancer 1 protein (AIB-1) CBP-interacting protein (pCIP) Class E basic helix-loop-helix protein 42 (bHLHe42) Receptor-associated coactivator 3 (RAC-3) Steroid receptor coactivator protein 3 (SRC-3) Thyroid hormone receptor activator molecule 1 (TRAM-1) NCOA3 AIB1 BHLHE42 RAC3 TRAM1
Observed MW:	157 kDa
Uniprot ID :	Q9Y6Q9

Application

Reactivity:	Human
Tested Application:	ELISA, WB, FC
Recommended dilution:WB:	1:500-1:1000
Image:	

Raji cells were subjected to SDS PAGE followed by western blot with FNab08218(SRC-3/AIB1 antibody) at dilution of 1:500

