

SRFBP1 antibody

Product Information

Catalog No.:	FNab08223
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in regulating transcriptional activation of cardiac genes during the aging process. May play a role in biosynthesis and/or processing of SLC2A4 in adipose cells(By similarity).

Immunogen information


Immunogen:	serum response factor binding protein 1
Synonyms:	Serum response factor-binding protein 1 SRF-dependent transcription regulation-associated protein p49/STRAP SRFBP1
Observed MW:	50 kDa
Uniprot ID :	Q8NEF9

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab08223(SRFBP1 antibody) at dilution of 1:100



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab08223(SRFBP1 Antibody) at dilution of 1:300