

SRSF3 antibody

Product Information

Catalog No.:	FNab08231
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	splicing factor, arginine/serine-rich 3
Synonyms:	Serine/arginine-rich splicing factor 3 Pre-mRNA-splicing factor SRP20 Splicing factor, arginine/serine-rich 3 SRSF3 SFRS3 SRP20
Observed MW:	22 kDa
Uniprot ID :	P84103

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse testis tissue slide using FNab08231(SRP20 Antibody) at dilution of 1:50



Various lysates were subjected to SDS PAGE followed by western blot with FNab08231(SRP20 antibody) at dilution of 1:1000