

## SSH3 antibody

### Product Information

Catalog No.:	FNab08248
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Protein phosphatase which may play a role in the regulation of actin filament dynamics. Can dephosphorylate and activate the actin binding/depolymerizing factor cofilin, which subsequently binds to actin filaments and stimulates their disassembly(By similarity).

### Immunogen information


Immunogen:	slingshot homolog 3(Drosophila)
Synonyms:	Protein phosphatase Slingshot homolog 3 SSH-like protein 3 (SSH-3L, hSSH-3L) SSH3 SSH3L
Observed MW:	90-95 kDa
Uniprot ID :	Q8TE77

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab08248(SSH3 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08248(SSH3 antibody) at dilution of 1:800