

SSR4 antibody

Product Information

Catalog No.:	FNab08249
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


TRAP proteins are part of a complex whose function is to bind calcium to the ER membrane and thereby regulate the retention of ER resident proteins.

Immunogen information


Immunogen:	signal sequence receptor, delta(translocon-associated protein delta)
Synonyms:	Translocon-associated protein subunit delta (TRAP-delta) Signal sequence receptor subunit delta (SSR-delta) SSR4 TRAPD
Observed MW:	19 kDa
Uniprot ID :	P51571

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab08249(SSR4 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08249(SSR4 antibody) at dilution of 1:500