

SSTR3 antibody

Product Information

Catalog No.:	FNab08254
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Receptor for somatostatin-14 and-28. This receptor is coupled via pertussis toxin sensitive G proteins to inhibition of adenylyl cyclase.

Immunogen information

Immunogen:	somatostatin receptor 3
Synonyms:	Somatostatin receptor type 3 (SS-3-R, SS3-R, SS3R, SST3) SSR-28 SSTR3
Observed MW:	None
Uniprot ID :	P32745

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab08254(SSTR3 antibody) at dilution of 1:50