

SSTR5 antibody

Product Information

Catalog No.:	FNab08255
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Receptor for somatostatin 28 and to a lesser extent for somatostatin-14. The activity of this receptor is mediated by G proteins which inhibit adenylyl cyclase. Increases cell growth inhibition activity of SSTR2 following heterodimerization.

Immunogen information

Immunogen:	somatostatin receptor 5
Synonyms:	Somatostatin receptor type 5 (SS-5-R, SS5-R, SS5R, SST5) SSTR5
Observed MW:	50 kDa, 65 kDa
Uniprot ID :	P35346

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human placenta tissue were subjected to SDS
PAGE followed by western blot with
FNab08255(SSTR5 antibody) at dilution of 1:1200