

SSX1 antibody

Product Information

Catalog No.:	FNab08257
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Could act as a modulator of transcription.

Immunogen information

Immunogen:	synovial sarcoma, X breakpoint 1
Synonyms:	Protein SSX1 Cancer/testis antigen 5.1 (CT5.1) Synovial sarcoma, X breakpoint 1 SSX1
Observed MW:	21-25 kDa
Uniprot ID :	Q16384

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab08257(SSX1 Antibody) at dilution of 1:300