

ST8SIA2-Specific antibody

Product Information

Catalog No.:	FNab08276
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ST8 alpha-N-acetyl-neuraminide alpha-2,8-sialyltransferase 2
Synonyms:	Alpha-2,8-sialyltransferase 8B Sialyltransferase 8B (SIAT8-B) Sialyltransferase St8Sia II (ST8SiaII) Sialyltransferase X (STX) ST8SIA2 SIAT8B STX
Observed MW:	42 kDa
Uniprot ID :	Q92186

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab08276(ST8SIA2-Specific antibody) at dilution of 1:50



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab08276(ST8SIA2-Specific antibody) at dilution of 1:400