

## STAC2 antibody

### Product Information

Catalog No.:	FNab08279
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	SH3 and cysteine rich domain 2
Synonyms:	SH3 and cysteine-rich domain-containing protein 2 24b2/STAC2 Src homology 3 and cysteine-rich domain-containing protein 2 STAC2
Observed MW:	42-45 kDa
Uniprot ID :	Q6ZMT1

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate hyperplasia slide using FNab08279(STAC2 Antibody) at dilution of 1:50



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab08279(STAC2 Antibody) at dilution of 1:300