

## STARD3NL antibody

### Product Information

Catalog No.:	FNab08290
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	STARD3 N-terminal like
Synonyms:	STARD3 N-terminal-like protein MLN64 N-terminal domain homolog STARD3NL MENTHO
Observed MW:	27 kDa
Uniprot ID :	O95772

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab08290(STARD3NL antibody) at dilution of 1:100



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab08290(STARD3NL antibody) at dilution of 1:1000