

## STAT4 antibody

### Product Information

Catalog No.:	FNab08302
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Carries out a dual function: signal transduction and activation of transcription. Involved in IL12 signaling.

### Immunogen information


Immunogen:	signal transducer and activator of transcription 4
Synonyms:	Signal transducer and activator of transcription 4 STAT4
Observed MW:	85 kDa
Uniprot ID :	Q14765

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:500
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab08302(STAT4 Antibody) at dilution of 1:200



Hela cells were subjected to SDS PAGE followed by western blot with FNab08302(STAT4 antibody) at dilution of 1:1000