

STAT6 antibody

Product Information

Catalog No.:	FNab08307
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Carries out a dual function: signal transduction and activation of transcription. Involved in IL4/interleukin-4-and IL3/interleukin-3-mediated signaling.

Immunogen information


Immunogen:	signal transducer and activator of transcription 6, interleukin-4 induced
Synonyms:	Signal transducer and activator of transcription 6 IL-4 Stat STAT6
Observed MW:	100-110 kDa
Uniprot ID :	P42226

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human spleen using FNab08307(STAT6 antibody) at dilution of 1:50



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab08307(STAT6 antibody) at dilution of 1:1500