

STAU2 antibody

Product Information

Catalog No.:	FNab08313
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


RNA-binding protein required for the microtubule-dependent transport of neuronal RNA from the cell body to the dendrite. As protein synthesis occurs within the dendrite, the localization of specific mRNAs to dendrites may be a prerequisite for neurite outgrowth and plasticity at sites distant from the cell body(By similarity).

Immunogen information

Immunogen:	stau2, RNA binding protein, homolog 2(Drosophila)
Synonyms:	Double-stranded RNA-binding protein Stau2 homolog 2 STAU2
Observed MW:	62-66 kDa
Uniprot ID :	Q9NUL3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab08313(STAU2 antibody) at dilution of 1:1000